

SC²P**Standard Calcium & Channel Screening**(this protocol must be used with the C²EP protocol A&B)

Notes: To ensure a fast, reproducible and fast screening of transfected/infected cells on Ca²⁺/channel signaling the following standard protocols should routinely be repeated at least on three different days or with at least three different cell patches.

Data Acquisition:

1. Make sure that you received all information necessary to complete the CCEP. No exceptions!
2. Carefully explain the experimental parameters. The time you win not do so is minimal but you save hours or days afterwards.


Data Analysis:

1. Perform data analysis immediately!
2. Analyse first, be critical to our hypothesis and try to extract as many as possible information out of the experiments (e.g. basal Ca²⁺, amplitude of release, area under the curve, time correlation, etc.). Never forget that subpopulations are possible particular in transfected cells (e.g. different amount of protein expression).

Standard Screening Protocols**Standard Protocol 1:****"Normal Ca²⁺ protocol"**

Puffers used: 2NaCa 0CaNa 0CaNa, EGTA

After 1 min experimental time to monitor basal Ca²⁺ concentration add the agonist to be tested (e.g. 10 μM histamine). Wait until a constant plateau or a clear and constant decline to basal level is shown. Do not stop too early - You might lose information!


**Data analysis:**

- Check basal Ca²⁺ levels (is it constant or changing?)
- Calculate raise in Ca²⁺ in response to the agonist. Don't forget to compare even the onset time and speed of elevation between given groups.
- Calculate Ca²⁺ plateau (if there is one) or the time constant for decline (τ = half time until the Ca²⁺ concentration reaches the resting level)
- Any oscillations or other events during onset, plateau or decline?

Standard Protocol 2:**"Ca²⁺ Oscillation protocol"**

Puffers used: 2NaCa 0CaNa 0CaNa, EGTA

After 1 min in HEPES+2.5 mM Ca_e switch to Ca²⁺-free solution and monitor basal Ca²⁺ over additional 1 min. Thereafter add the agonist (e.g. 10 μM histamine) for 2 to 3 min. Wait until oscillation stops or Ca²⁺ returns to basal. Do not stop too early - You might lose information!

**Data analysis:**

- Check basal Ca²⁺ levels (is it constant or changing?)
- Calculate 1st raise in Ca²⁺ in response to the agonist. Don't forget to compare even the onset time and speed of elevation between given groups.
- Calculate Ca²⁺ threshold for oscillation (if there is one), the amplitude and the frequency.
- Check the kinetics of oscillation by duration of Ca²⁺ elevation, τ , etc.


Standard Protocol 3:**"Ca²⁺ Release protocol"**

Puffers used: ✓ 2NaCa

○ 0CaNa

✓ 0CaNa, EGTA

After 1 min in HEPES+2.5 mM Ca_e switch to Ca²⁺-free solution and monitor basal Ca²⁺ over additional 1 min. Thereafter add the agonist (e.g. 10 μM histamine) for 1 to 2 min. Wait until Ca²⁺ returns to basal level (app. 3 to 4min) and add HEPES + 2.5 mM Ca_e (app. time point 6 min) Do not stop too early! Wait until after addition of Ca_e a plateau has been reached or Ca returns to baseline!

**Data analysis:**

- Check basal Ca²⁺ levels (is it constant or changing?)
- Calculate Ca²⁺ transient in response to the agonist. Don't forget to compare even the onset time, area under the curve and speed of elevation/decline between given groups.
- Calculate the raise in Ca²⁺ in response to the addition of Ca_e. Check the kinetics of this increase and whether plateau or decline.


Standard Protocol 4:**"Ca²⁺ Leak protocol"**

Puffers used: ✓ 2NaCa

○ 0CaNa

✓ 0CaNa, EGTA

After 1 min in HEPES+2.5 mM Ca_e switch to Ca²⁺-free solution and monitor basal Ca²⁺ over additional 4 to 5 min. Add HEPES + 2.5 mM Ca_e (app. time point 5 to 7 min). Do not stop too early! You might lose information!

**Data analysis:**

- Check basal Ca²⁺ levels (is it constant or changing?)
- Calculate raise in Ca²⁺ in response to the addition of Ca_e. Check the kinetics of this increase and whether plateau or decline.

Data process and handling:

1. Store your original printouts and the protocols together with the analysis of the data.
2. Provide WFG a Prizm printout and the file with your analysis via e-mail.
3. Data analysis must be provided with given mean ± sem (and n!).
4. Provide statistical evaluation using the Prizm program (unpaired t-test). Level of significance should be hold at p<0.05.
5. Stay in contact with other team members on their effort. Talk and discuss your findings with all the people in the team and tell them your findings. Do not hide but show your data. Just discussion makes us using forward.